

RECEIVED
CENTRAL TAX CENTER
NOV 26 2008

In the Specification:

Page 26, rewrite lines 3 to 13 as follows:

The reaction conditions include 0.2 μ M of Fwd1 for the products of 3' specific cDNA and 0.2 μ M of Rev1 for the products of 5' specific cDNA, 200 μ M dNTPs, 40 mM Tricine-KOH (pH 8.7), 15 mM KOAc, 3.5 mM Mg(Oac)₂, 3.75 g/ml BSA, 0.005% Tween-20 (polyacetate), 0.005% Nonidet-P40, and 0.5 U Taq DNA polymerase in a final volume of 50 μ l. The PCR reactions are carried out in a Perkin-Elmer 9700 thermocycler using the following thermal cycle parameters: 5 cycles comprising a denaturation at 94°C for 5 seconds, a hybridization of the primers at 72°C, 5 cycles comprising a denaturation at 94°C for 5 seconds, a hybridization of the primers at 70°C for 10 seconds, and an extension of polymerization at 72°C for 3 minutes and finally 25 cycles comprising a denaturation at 94°C for 5 seconds, a hybridization of the primers at 68°C for 10 seconds, and a polymerase extension at 72°C for 3 minutes.

Page 38, rewrite lines 1 to 11 as follows:

The part of the cDNA encoding for heterocarpine is inserted at the BamH1/XhoI sites of the pQE-TriSystem (Qiagen) expression vector. The pQE-TriSystem vector contains the activating sequences of the cytomegalovirus (CMV) fused to chicken beta-actin promoter allowing a very significant heterologous expression. Human embryo kidney (HEK-293) cells are cultured in DMEM medium (Dulbecco's Modified Eagle's Medium) containing

100 U/ml of penicillin and 100 µg/ml of strptomycin sulphate, complemented with 10% foetal calf serum. The cells are sub-cultured 24 hours before the transfection protocol allowing normal metabolism of the cells and better transfection efficiency. The transfection of 1 µg of pQE-TruSystem containing the cDNA encoding for heterocarpine was carried out using Effectene® Transfection reagent according to the manufacturer's (Qiagen) recommendations.